Med-Fit 1

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Skin Preparation

Before Application of Self-Adhesive Electrodes

For best results please wash the area to be treated with a solution of warm salt water (mix one tablespoon of salt in warm water) and rub the area with the salt solution.

The procedure will maximise the effect of the Tens when applied to areas of the body, where the skin may be particularly dry.

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INDEX

Chapter	Contents	Page
1.	General Description	1
2.	Warnings & Precautions	2
3.	Control Layout.....	3
4.	Technical Specifications	4
5.	Accessories	5
6.	Graphic Symbols	6
7.	Parameter Controls	6
8.	Attachment of Electrodes Lead Wires	9
9.	Preparation for use	10
10.	Application of Re-usable Self Adhesive Electrodes	10
11.	Adjusting The Controls	11
12.	Battery Information	14
13.	Maintenance Transportation and Storage	16
14.	Safety-technical Controls.....	16
15.	Conformity To Safety Standards	17
16.	Warranty	17
	Manufacturer	17
	Representative In The EU	17

Chapter 1 : GENERAL DESCRIPTION

The MED-FIT 1 is a battery operated pulse generator that sends electrical impulses to the body by means of electrodes attached to the skin.

The device has two independently controllable channels. A pair of electrodes can be connected to each channel.

The electronics of the MED-FIT 1 create electrical impulses whose Intensity, duration, number per second and modulation may be altered with the controls or switches. Dial controls are very easy to use and

the slide cover prevents accidental changes to the settings.

Chapter 2: WARNINGS & PRECAUTIONS

PLEASE NOTE:

It is imperative that patients read and understand the warnings and precautions before using the unit. Do not allow your machine or electrodes to be used by anyone else, as they are designed for single patients use only. It is recommended that proper medical advice on the use of TENS is sought from a Qualified Practitioner (Physiotherapist, Doctor or Nurse) prior to use, in order to ensure safe and effective treatment. If you are taking any medication please carry on as normal but seek advice from your clinician before using the device.

WARNING! PATIENTS WITH PACEMAKER MAY NOT BE TREATED WITH TENS

- Do not use during pregnancy except during labour
- Do Not place electrodes over the Carotid Sinus
- Do Not use on broken or damaged skin
- Do Not place electrodes close to the eyes or in the mouth.
- Do Not use TENS whilst driving or operating machinery.

Tens is unsuitable and should not be used in the following situations.

- Persons suffering from conditions where the circulation is impaired.
- Epilepsy, Heart Condition or any form of Malignancy.
- Patients with poor skin sensation and non-compliant patients who are emotionally disturbed or have dementia.
- Over metal implants or in conjunction with sleep apnea or heart monitors.

You should be aware that TENS units provide symptomatic relief only and are not considered curative.

The degree of pain relief declines in some cases with time, TENS may have less effect after extended use. To reduce the risk of development of tolerance, the following measures may be taken:

- Frequent follow-ups.
- Teach the patient to use burst, modulated, and conventional stimulation.
- Vary electrode placement.
- Vary the frequency


If none of the above are successful discontinue TENS for up to 2 weeks and then resume stimulation.

As with pain-killing medicines, TENS provides temporary pain relief, which often lasts up to four hours after treatment. Treatment usually lasts 30-60 minutes, 2-4 times daily.

The time taken to achieve pain relief varies from an immediate effect up to about an hour. Some patients prefer continuous stimulation while others prefer intermittent stimulation with different intervals.

Remove the electrodes and discontinue stimulation if you experience skin irritation or discomfort.

Chapter 3 : CONTROL LAYOUT

Chapter 4 : TECHNICAL SPECIFICATIONS

The technical specification details of MED-FIT 1 are as follows

	MECHANISM	TECHICAL DESCRIPTION
01	Channel	Dual, isolated between channels
02	Pulse Amplitude	Adjustable 0-80mA, Max output 80mA peak to peak (15.8mA rms) into 500ohm load each channel.
03	Output Voltage	Adjustable 0-40V, Max output 40V peak to peak (7.9V rms) into 500ohm load each channel.
04	Wave Form	Asymmetrical Bi-Phasic Square Pulse
05	Power supply	One 9 Volt Battery, type 6F22
06	Size	95(H) x 65(W) x 23.5(T) mm
07	Weight	115 grams (battery included)
08	Pulse Rate	Adjustable, from 2 Hz to 150 Hz
09	Pulse Width	Adjustable, from 30µs to 260µs
10	Burst Mode (B Mode)	Bursts occur twice per second. The pulse width is adjustable. The frequency is fixed at 100Hz
11	Modulation	Pulse rate is automatically varied in a cyclic pattern over a nominal 10 second interval (in max 150Hz). Pulse rate decreases linearly, over a period of 4 seconds, from the control setting value to a value which is 40% less. The lower pulse rate will continue for 1 second and then the pulse rate increases linearly over a period of 4 seconds to the original value. This will continue for 1 second and the cycle is then repeated
12	Operating Condition	Temperature : 0°~40°C Relative Humidity : 30%~75% Atmosphere Pressure : 700Hpa~1060Hpa
13	Remark	There may be up to a +/- 5% tolerance of pulse width and pulse rate and +/- 20% tolerance of amplitude & voltage.


Chapter 5 : ACCESSORIES

Each set MED-FIT 1 are completed with standard accessories and standard label as given below

I. Accessories

REF. NO.	DESCRIPTION	Q'TY
1. SA 10	50 X 50 mm Adhesive Electrodes	4 pieces
2.	Electrodes Leads	2 pieces
3.	9 V Battery, type 6F22	1 piece
4.	Instruction Manual	1 piece

II. LABEL

The label attached to the back of the unit contains important information this model, serial number, supply voltage, the name of manufacturer, CE number and classification. Please do not remove.

Chapter 6 : GRAPHIC SYMBOLS

1.  Note Operating Instructions
2.  Degree of Electrical Protection BF
3.  Do not insert the plug into AC power supply socket.
4.  Direct Current (DC power source)
5.  Consult Instructions for use
6.  Manufacturer
7.  Serial Number

Chapter 7 : PARAMETER CONTROLS

PULSE DURATION

Wider pulse duration settings will deliver stronger stimulation for any given intensity setting. As mentioned in the Controls section, by using a combination of intensity and pulse duration different groups of nerve fibres can be stimulated.

The wider pulse duration is needed to recruit motor fibres, whereas the narrow pulse duration is used on the sensory fibres.

The choice of which pulse duration to use is partially dependent upon the Treatment Mode and Protocol selected (refer to the appropriate section).

PULSE RATE

The Pulse Rate (frequency, pulses per second or Hertz) is thought to determine whether the spinal gating mechanism or the endorphin release mechanism is activated.

As a rough guide the higher frequencies, those above 70Hz, will activate the spinal gating mechanism. The patient will perceive a "pins and needles" sensation.

Lower frequencies, those between 2-16Hz, will activate the endorphin release mechanism. The patient will perceive a "tapping" sensation.

It should be noted that individual patients vary in their response to TENS and the practitioner should be prepared to experiment with different settings in order to obtain the most effective response.

TREATMENT MODE

Normal or Conventional TENS offers the practitioners complete control over all the various treatment parameters of the instrument.

Burst Mode is analogous to the Low Rate TENS technique except the low frequency individual pulses are replaced by individual "bursts" of 7-10 individual pulses. It is thus a combination of Conventional TENS and Low Rate TENS. In Burst Mode, the treatment frequency is fixed by the instrument and is not adjustable with the Frequency Rate control.

Modulated Mode attempts to prevent nerve accommodation by continuously cycling the treatment intensity. When using Modulated Mode increase the intensity only when the unit is at the maximum intensity of the modulation cycle. If the intensity is increased during a low intensity period of the cycle the resulting stimulation may cause discomfort.

INTENSITY

Each patient responds differently to different levels of intensity, due to varying degrees of tissue resistance, innervation, skin thickness, etc. Intensity instructions are therefore limited to the following settings: Perception - The intensity is increased so that the patient can feel the stimulation, but there is no muscular contraction.

Slight Contraction - Intensity is increased to a barely visible muscular contraction that is not strong enough to move a joint. When using low pulse rate settings, this will show as individual twitches. At higher pulse rates there will simply be increased muscle tension. Strong muscular contraction is typically not used in TENS therapy. However, muscular contraction may be useful if the pain involves a cramped or spastic muscle. The TENS can be used as a traditional muscle stimulator in the circumstances to quickly break the spasm. Use a rapid pulse rate, wide pulse duration and set the intensity to visible contraction (still within patient tolerance). Twenty or thirty minutes of such a tetanized muscular contraction will generally break the spasm. In all cases, if the patient complains that the stimulation is

uncomfortable, reduce intensity and/or cease stimulation.

TIME DURATION

The onset of pain relief should occur shortly after the intensity setting has been determined. However, in some cases, pain relief may take as long as 30 minutes to achieve especially when using low-frequency TENS. TENS units are typically operated for long periods of time, with a minimum of 20 - 30 minutes and in some post-operative protocols, as long as 36 hours.

In general, pain relief will diminish within 30 minutes of the cessation of stimulation. Pain relief obtained through low-frequency TENS may last longer.

Chapter 8 : ATTACHMENT OF ELECTRODE

LEAD WIRES

The wires provided with the system insert into the jack sockets located on top of the device. Holding the insulated portion of the connector, push the plug into one of the jacks (see drawing); one or two sets of wires may be used.

After connecting the wires to the stimulator, attach each wire to an electrode. Take care when you attach and detach the wires. Jerking the wire instead of holding the insulated connector body may cause wire breakage.

CAUTION

Do not insert the plug of the patient lead wire into the AC power supply socket.

Chapter 9 : PREPARATION FOR USE

It is recommended that for persons with sensitive skin, the area is prepared prior to the application of the self adhesive electrodes by cleaning using soap and water or alcohol wipes and then drying thoroughly. It is quite normal to see a reddening of the skin underneath the electrodes following treatment but this usually disappears in less than 2 hours.

If the reddening persists we recommend PALS (BLUE) electrodes which are specifically designed for sensitive skin.

Chapter 10 : APPLICATION OF RE-USABLE SELF ADHESIVE ELECTRODES

Application

1. Clean and dry the skin at the thoroughly with soap and water prior to application of electrodes.
2. Insert the lead wire into the pin connector on the pre-wired electrodes.
3. Remove the electrodes from the protective liner and apply them firmly to the treatment site.

Removal


1. Peel the electrodes from the skin; do not pull on the lead wires because it may damage the electrodes.
2. Place the electrodes on the liner and remove the lead wire by gently twisting and pulling.

Care and Storage

1. Between uses, store the electrodes on their liner in the resealable bag in a cool dry place.
2. The adhesive properties of the electrodes may be improved by applying a few drops of cold water and allowing the electrode surface to dry. Do NOT soak the electrodes.

Important

1. Do not apply to broken skin.
2. The electrodes should be discarded when they no longer adhere.
3. The electrodes are intended for single patient use only.
4. If irritation occurs, discontinue use and consult your clinician.
5. Read the instruction for use of self-adhesive electrodes before application.


Chapter 11 : ADJUSTING THE CONTROLS

1. Slide Cover :

A sliding panel covers the controls for Pulse Width, Pulse Rate, Mode Selector and Modulation Selector. Take advice from your clinician regarding the settings and then leave the cover in place.


2. Display LED

The LED corresponding to the appropriate channel will light whenever that channel is activated. At low frequencies each individual pulse may be observed but at frequencies above 30 Hz the LED will appear to be constantly illuminated.

3.On/Off Switch and Intensity Control :

The unit is switched on by turning the appropriate channel control clockwise. The LED will illuminate as described above. The output intensity increases as the control is turned further clockwise. To reduce the output intensity and/or switch off the unit the control is turned anti-clockwise.


4.Lead Connector

Ensure that the unit is switched off. Connect the electrodes to the lead by inserting the plugs into the electrode connectors. Plug the lead into the unit and press the electrodes firmly onto the skin.


5. Mode Control

Expose the controls by sliding front cover down. This switch has 3 positions: B for Burst stimulation, N for Constant stimulation, and M for modulated stimulation. Push the Mode Selector to the desired position.


6. Pulse Rate Control:

This control determines the frequency (pulses per second). This dial controls the frequency to both channels. Unless otherwise instructed set the pulse rate to between 70 - 120 Hz.

7. Pulse Width Control:


This control adjusts the length of each electrical impulse. By increasing the pulse width the strength and sensation of the stimulus will be increased. By reducing the pulse width the strength and sensation will be reduced. Unless otherwise instructed set the control to between 70 - 120 μ s.

8.Check/Replace the Battery:

Over time, in order to ensure the functional safety of TENS, changing the battery is necessary.

1. Make sure that both intensity controls are in the off position
2. Slide the battery compartment cover and remove.
3. Remove the battery from the compartment.
4. Insert the battery into the compartment. Note the polarity indicated on the battery and in the compartment.
5. Replace the battery compartment cover and slide to close

Chapter 12 : BATTERY INFORMATION

PRECAUTIONS

1. Remove battery if equipment is not likely to be used for some time.
2. Please recycle the used battery in accordance with domestic regulation.
3. Do not throw the used battery into fire.

If you use rechargeable batteries, please follow the instructions.

RECHARGEABLE BATTERIES(NOT INCLUDED)

Prior to the use of a new unit, the rechargeable battery should be charged according to the battery manufacturer's instructions. Before using the battery charger read all the instructions and cautionary markings on the battery and in this instruction manual.

After being stored for 60 days or more, the batteries may lose their charge. After long periods of storage, batteries should be charged prior to use.

BATTERY CHARGING

- (1) Plug the charger into any working 110 or 220/240v mains electrical outlet. The use of any attachment not supplied with the charger may result in the risk of fire, electric shock, or injury to persons.
- (2) Follow the battery manufacturer's instructions for charging time.
- (3) After the battery manufacturer's recommended charging time has been completed, unplug the charger and remove the battery.
- (4) Batteries should always be stored in a fully charged state.

To ensure optimum battery performance, follow these guidelines:

- (a) Although overcharging the batteries for up to 24 hours will not damage them, repeated overcharging may decrease useful battery life.
- (b) Always store batteries in their charged condition. After a battery has been discharged, recharge it as soon as possible. If the battery is stored more than 60 days, it may need to be recharged.
- (c) Do not short the terminals of the battery. This will cause the battery to get hot and can cause permanent damage. Avoid storing the batteries in your pocket or purse where the terminals may accidentally come into contact with coins, keys or any metal objects.
- (d) **WARNINGS:**
 1. Do not attempt to charge any other types of batteries in your charger, other than the nickel-cadmium rechargeable batteries. Other types of batteries may leak or burst.
 2. Do not incinerate the rechargeable battery as it may explode!

Chapter 13 : MAINTENANCE, TRANSPORTATION AND STORAGE OF TENS DEVICE

1. Non-flammable cleaning solution is suitable for cleaning the device.
Note: Do not smoke or work with open lights (for example, candles, etc.) when working with flammable liquids.
2. Stains and spots can be removed with a cleaning agent.
3. Do not submerge the device in liquids or expose it to large amounts of water.
4. Return the device to the carrying box with sponge foam to ensure that the unit is well-protected before transportation.
5. If the device is not to be used for a long period of time, remove the batteries from the battery compartment (acid may leak from used batteries and damage the device). Put the device and accessories in carrying box and keep it in cool dry place.
6. The packed EMS device should be stored and transported under the temperature range of -20°C ~ + 60°C, relative humidity 20% ~ 95%, Atmosphere pressure 500hPa ~ 1060hPa.

Chapter 14 : SAFETY-TECHNICAL CONTROLS

For safety reasons, check your MED-FIT 1 each week based on the following checklist.

1. Check the device for external damage.
 - deformation of the housing.
 - damaged or defective output sockets.
2. Check the device for defective operating elements.
 - legibility of inscriptions and labels.
 - make sure the inscriptions and labels are not distorted.
3. Check Led
 - led must be illuminated when switched on.
4. Check the usability of accessories.
 - patient cable undamaged.
 - electrodes undamaged.

Please consult your distributor if there are any problems with device and accessories.

Chapter 15 : CONFORMITY TO SAFETY STANDARDS

The MED-FIT 1 devices are in compliance with the EN 60601-1-2:2001 and EN 60601-1:1990+A1:1993+A2:1995+A13:1996 safety standards.

Chapter 16 : WARRANTY

All MED-FIT 1 models carry a warranty of one year from the date of delivery. The warranty applies to the stimulator only and covers both parts and labour relating thereto.

The warranty does not apply to damage resulting from failure to follow the operating instructions, accidents, abuse, alteration or disassembly by unauthorized personnel.

Manufacturer:

Everyway Medical Instruments Co., Ltd.
3F., No.5, Ln. 155, Sec. 3, Beishen Rd.,
Shenkeng Dist., New Taipei City 22203,
Taiwan. (R.O.C.)

Representative in the EU:

REHAB EUROPA SL
SANT GERVASI DE CASSOLES, 96 3º 4^a
08022 BARCELONA, SPAIN.

INFORMATION FOR DISTRIBUTOR:

Please contact the above mentioned manufacturer for technical support and documentation when necessary.

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Copyright©2011 by Everyway Medical Instruments Co., Ltd.; Edition: V1.4
Med-fit Ltd