Safe-Trac Overdoor Cervical Traction

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
General and Assembly Instructions

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Instructions for home use

Consult a health professional before using
for the first time

- Fit the rubber pad into the door bracket this protects the door and also gives a good fit
- Attach bracket to door
- Attach the traction unit to the bracket as shown Fig A
- Place a chair against the door ready for use

- Fit the harness to the back of the head and under the chin
- Fasten the side straps on the sides of the chin ensuring a tight fit
- Pull the cord and wait for a click. Assess whether there is any tension or traction.
- Pull the cord again and wait for the next click. Pause again and then repeat.
- Once tension or traction is detected begin to count the clicks as each click from this point on denotes 1lb/0.45kg of traction force.

- When time is up release the side straps from the chin strap.

DO NOT REACH UP TO RELEASE THE TAB ON THE RATCHET.

Recommended use: 5 – 10 clicks (equates to 5 – 1-lbs/2.26 – 4.53kg of force) for 10 minutes morning and evening until pain relief is felt. Continue to use as required

WARNING: Consult a health professional before using for the first time

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Med-Fit UK Ltd